


PAJUNK®



UniPlex®

*Stimulating Single Shot
Nerve Block Needles*



MADE IN GERMANY 

The Stimulating Needle with NanoLine® Coating

The UniPlex NanoLine needles have been developed for single shot nerve block applications, allowing for stimulation guidance only. With the introduction of the NanoLine thin layer coating technology, PAJUNK® became the innovative technology leader in stimulation guided nerve block needles. The new type of coating offers improved gliding properties and precise stimulation while reducing needle thickness.

→ All coated Pajunk nerve block needles use NanoLine technology offering superior gliding properties



AVAILABLE OPTIONS

UniPlex II NRFit®

Stimulating single shot nerve block needles with non-removable injection tubing, removable stimulation cable, colour-coded needle hub and NanoLine Coating

More information at: pajunk.com



NanoLine Coating

Very thin polymer film, smooth surface,
accurate inner and outer diameter

- Excellent puncture and gliding properties through smooth surface
- Improves visibility under ultrasound¹
- Precise stimulation through the non-insulated needle tip

Depth Markings

- Easy to read and identify

Echogenic Needle Tip

Facet tip with two inclination angles

- Improves needle tip visibility under ultrasound

Facet tip



Also available in NRFit®

UniPlex needle
with facet tip

Size	UniPlex Item no.	UniPlex II NRFit® Item no.	PU
25G x 35 mm (1 1/3")	001156-80	001246-80	10
25G x 50 mm (2")	001156-81	001246-81	10
22G x 50 mm (2")	001156-74	001246-74	10
22G x 80 mm (3 1/5")	001156-71	001246-71	10
22G x 100 mm (4")	001156-84		10
22G x 120 mm (4 3/4")	001156-82	001246-82	10
21G x 100 mm (4")	001156-77	001246-77	10
20G x 120 mm (4 3/4")	001156-72	001246-72	10
20G x 150 mm (6")	001156-76	001246-76	10

UniPlex needle
with facet S tip

Size	UniPlex Item no.	UniPlex II NRFit Item no.	PU
24G x 50 mm (2")	001157-81	001247-85	10
22G x 50 mm (2")	001157-74	001247-74	10
22G x 80 mm (3 1/5")	001157-71		10
21G x 100 mm (4")	001157-77		10

UniPlex needle
with SPROTTE® tip

Size	UniPlex Item no.	UniPlex II NRFit Item no.	PU
22G x 50 mm (2")	001156-31G	001246-31G	10
22G x 70 mm (2 3/4")	001156-31H	001246-31H	10