


Manometer

Manometer for CSF Pressure Measurement

As the manufacturer of the world's first atraumatic spinal needle, Pajunk takes the diagnostic spinal tap as an integral part of its core competences. With the SPROTTE® Lumbar needle and our Manometer, Pajunk provides the complete equipment for diagnostic CSF sampling and pressure measurement. The Pajunk spinal manometer allows for measurements up to 54 cm H₂O (with the additional extension tubing) and is available in Luer and NRFit®.



- A solid, transparent tubing allows for easy measurements up to 34cm H₂O
- Easy and secure assembly with Luer or NRFit connection
- Extension tube available separately allowing for measurements up to 54cm H₂O

Manometer	Item description	Item no.	NRFit	PU
			Item no.	
	3-34 cm H ₂ O	001151-38G	001163-38F	10
	Extension Tubing 20cm H ₂ O	001152-38G	001152-38G	10

Injection Tube NRFit male/ NRFit female	Item description	Item no.	NRFit	PU
			Item no.	
		20cm (8")		15461020
	60cm (24")		001163-37R	10

Injection Tube LLA male / LLA female	Item description	Item no.	NRFit	PU
			Item no.	
	60cm (24")	001151-37R		10

OPTIDYNAMIC® Spinal Manometer

For diagnostic CSF pressure measurement

The MediPlast OPTIDYNAMIC Spinal Manometer is characterized by its intuitive and simple handling. The very precise CSF pressure measurement is performed by a thin and flexible tube which enables the use of the manometer in different patient positionings.

The combination of a PAJUNK SPROTTE® cannula for lumbar puncture and the Optidynamic Spinal Manometer thus enables secure and patient-friendly diagnostics.



- Thin, calibrated tube for precise measurement performance at low liquor losses
- Flexible connection tubing for use in various storage positions
- Simple and intuitive handling

OPTIDYNAMIC Spinal Manometer

Item no.	NRFit® Item no.	PU
6061650008		25
	6061650009	20

PAJUNK® GmbH Medizintechnologie
Global Headquarters
Karl-Hall-Straße 1
78187 Geisingen, Germany
Tel.: +49 7704 9291-0
Fax: +49 7704 9291-600
info@pajunk.com
pajunk.com

PAJUNK® Medical Produkte GmbH
Central Europe
Karl-Hall-Straße 1
78187 Geisingen, Germany
Tel.: +49 7704 8008-0
Fax: +49 7704 8008-150
info@pajunk-medipro.com
pajunk.de & pajunk.eu

PAJUNK® Swiss AG
Switzerland
Neuhofstrasse 5A
6340 Baar, Switzerland
Tel.: +41 41 760 96 27
info@pajunk.ch
pajunk.ch

PAJUNK® Medical System L.P.
United States of America
4575 Marconi Drive
Alpharetta, GA 30005, USA
Tel.: +1 770.493.6832
Fax: +1 678.514.3388
info@pajunk-usa.com
pajunkusa.com

PAJUNK® UK Medical Products Ltd
United Kingdom
Unit D1 The Waterfront
Goldcrest Way, Newburn Riverside
Newcastle upon Tyne, NE15 8NY, UK
Tel.: +44 191 264 7333
info@pajunk.co.uk
pajunk.co.uk

