


EpiLong II

*Catheters and Needles
for Epidural Anaesthesia*




Tuohy II Epidural Injection Needles

Following the original Tuohy tip design and made of high-quality polished steel, the PAJUNK® Tuohy needles are perfectly suitable for Epidural Anaesthesia.



Integrated Ergonomic Wings

→ Optimum handling for different gripping techniques



Colour-Coded Needle Hub

→ Easy identification of needle size



Highest Processing Quality

- High quality stainless steel needle for increased stability
- Consistent feel and reliable performance




Depth Markings

- Easy to read and identify



Tuohy tip



Precision-Ground Plastic Stylet

- Flush ends of plastic stylet and needle reduce the risk of tissue displacement
- Increased needle stability



Also available in **NRFit**

Epidural Anaesthesia with *EpiLong II & EpiLong II Soft*

PAJUNK® Epidural sets and trays come equipped with the Pajunk Tuohy needle and a selection of EpiLong II or EpiLong II Soft catheters as well as accessories like the clamping adapter, bacteria filters, and LOR syringes.




VARYING DEPTH MARKINGS BETWEEN LUER AND NRFIT®

- ▶ Due to the extended NRFit hub length markings on NRFit catheters are slightly different from the Luer versions.

Depth Markings


Graduations on the first 30 cm

- Easy determination of catheter positioning
- Easy to read and identify



Marking for Clamping Adapter

- Simplifies connection of catheter and adapter





Catheter with Closed Tip

→ With three lateral openings



Catheter with Closed Tip

→ With six lateral openings



Stainless Steel Helical Coil

→ Improved catheter flexibility and stability resulting in reduced chance of kinking¹

→ Radiopaque



Soft Tip

→ Minimizes the risk of accidental injury to the blood vessels and dura^{2,3}



Also available in **NRFit**

1. Toledano, R. D. et al. Epidural Catheter Design - History, Innovations, and Clinical Implications, Anesthesiology, V 121, No 1, 2014
 2. Banwell B. R. et al. Decreased incidence of complications in parturients with the Arrow (FlexTip Plus™) epidural catheter, CAN J ANAESTH 1998; 45:4; pp 370-72
 3. Juneja M. Incidence of Epidural Vein Cannulation in Parturients with Three Different Epidural Catheters, Regional Anesthesia Vol. 21 No. 2S March-April Supplement 1996

Tuohy Needles

Tuohy II
with plastic stylet and
integrated wings

Size Needle	Item no.	NRFit® Item no.	PU
19G x 90 mm (3 1/2")	1170-6I091		10
18G x 80 mm (3 1/4")	1170-6K081	1176-6K081	10
18G x 90 mm (3 1/2")	1170-6K091	1176-6K091	10
18G x 103 mm (4")	1170-6K103	1176-6K103	10
18G x 120 mm (4 3/4")	1170-6K120	1176-6K120	10
17G x 80 mm (3 1/4")	1170-6M081	1176-6M081	10
17G x 90 mm (3 1/2")	1170-6M091	1176-6M091	10
17G x 103 mm (4")	1170-6M103	1176-6M103	10
17G x 120 mm (4 3/4")	1170-6M120	1176-6M120	10

Tuohy
with metal stylet and
detachable wings

Size Needle	Item no.	NRFit® Item no.	PU
20G x 90 mm (3 1/2")	1150-4G090		25
19G x 50 mm (2")	1150-4I050		25
18G x 150 mm (6")	1150-4K150		10
17G x 150 mm (6")	1150-4M150		10

Caudal Needles

Caudal
with Crawford tip

Size Needle	Item no.	NRFit® Item no.	PU
25G x 25 mm (1")		1163-12C025	25
23G x 25 mm (1")	1149-12D025		25
22G x 40 mm (1 5/8")	1149-12E040	1163-12E040	25
22G x 60 mm (2 3/8")	1149-12E060		25

EpiLong Catheters

EpiLong

Size Needle	Closed/ Open Tip	Lateral Openings	Item no.	NRFit Item no.	PU	Steel mandrin	Stainless steel helical coil	Clamping Adapter
20G x 90 cm (35")	closed	3	001151-37C		10	•	•	•
20G x 90 cm (35")	closed	3	021151-37C	021163-37C	10			•
19G x 90 cm (35")	closed	3	001151-37A		10	•		•
19G x 90 cm (35")	closed	3		021163-37A	10			•

EpiLong Soft
with helical coil

Size Needle	Closed/ Open Tip	Lateral Openings	Item no.	NRFit Item no.	PU	Steel mandrin	Stainless steel helical coil	Clamping Adapter
20G x 90 cm (35")	closed	6	041151-37R	041163-37R	10		•	•
19G x 90 cm (35")	closed	6	041151-37P		10		•	•

EpiLong II Catheter Sets

	Size Needle	Size Catheter 3 lateral openings	Item no.	NRFit Item no.	PU	LOR	Clamping Adapter	Filter	FixoLong	FixoCath	NRFit Syringe & Filter needle	
EpiLong II	18G x 80 mm (3 1/4")	20G	033115301-51		10	Lock	•	•				
			053115301-51		10	Slip	•	•	•			
			033116301-51		10	Slip	•	•				
			133116301-51		10	Slip	•	•			•	
	18G x 90 mm (3 1/2")	20G	033115201-51		10	Lock	•	•				
			033125201-51		10	Lock	•	•		•		
			053115201-51		10	Slip	•	•	•			
			033116201-51		10	Slip	•	•			•	
	17G x 80 mm (3 1/4")	19G	033115301-48		10	Lock	•	•				
			033116301-48		10	Slip	•	•			•	
				133116301-48		10	Slip	•	•			•
	EpiLong Short Length	19G x 50 mm (2")	23G	0131153-50	0131163-50	10	Slip	•	•			
EpiLong II Soft	18G x 90 mm (3 1/2")	20G	053115701-51	053116701-51	10	Slip	•	•	•			
	17G x 90 mm (3 1/2")	19G	053115701-48	053116701-48	10	Slip	•	•	•			
EpiLong II VPC	18G x 90 mm (3 1/2")	20G	034115201-51		10	•	•	•				

Epidural Accessories

Accessories

Product	Description	Item no.	NRFit® Item no.	PU
Epilong VPC	10ml, Lock	0001151-49		10
LOR Syringe	Lock	0001151-66		10
LOR Syringe	Slip		0001166-65	10
Filter		001151-62	001163-62	10
FixoLong		001163-60	001163-60	10
FixoLong + Filter		001151-61	001163-61	10
Clamping Adapter	yellow, < 20G	001151-37V	001163-37V	10
Clamping Adapter	white, ≥ 20G	001151-38C	001163-38C	10
FixoCath			001151-37Z	10

NRFit® is a registered trademark of CEDSA used with their permission.



Scan now!

PAJUNK® GmbH Medizintechnologie
Global Headquarters
Karl-Hall-Straße 1
78187 Geisingen, Germany
Tel.: +49 7704 9291-0
Fax: +49 7704 9291-600
info@pajunk.com
pajunk.com

PAJUNK® Medical Produkte GmbH
Central Europe
Karl-Hall-Straße 1
78187 Geisingen, Germany
Tel.: +49 7704 8008-0
Fax: +49 7704 8008-150
info@pajunk-medipro.com
pajunk.de & pajunk.eu

PAJUNK® Swiss AG
Switzerland
Neuhofstrasse 5A
6340 Baar, Switzerland
Tel.: +41 417 609 627
info@pajunk.ch
pajunk.ch

PAJUNK® Medical System L.P.
United States of America
4575 Marconi Drive
Alpharetta, GA 30005, USA
Tel.: +1 888 972 5865
Fax: +1 678 514 3388
info@pajunk-usa.com
pajunkusa.com

PAJUNK® UK Medical Products Ltd
United Kingdom
Unit D1 The Waterfront
Goldcrest Way, Newburn Riverside
Newcastle upon Tyne, NE15 8NY, UK
Tel.: +44 191 264 7333
info@pajunk.co.uk
pajunk.co.uk