


## SonoBlock II

*Echogenic, Non-Stimulating  
Single Shot Plane Block Needle*

SonoBlock needles were developed by PAJUNK® for single shot nerve block applications without stimulating, allowing for placement under ultrasound guidance only. The needles also includes the patented Cornerstone Reflectors for optimized needle echogenicity.



<sup>1</sup> Edgcombe H., Hocking G. (2010 March–April) Sonographic identification of needle tip by specialists and novices: a blinded comparison of 5 regional block needles in fresh human cadavers. Reg. Anesth. Pain Med., 35(2): 207–211.

# SonoBlock II

Echogenic, Non-Stimulating  
Single Shot Plane Block Needle

SonoBlock II	Size	Item no.	NRFit® Item no.	PU
Facet tip	22G x 50mm (2")	001280-74	001260-74	10
	22G x 80mm (3 1/5")	001280-71	001260-71	10
	21G x 100mm (4")	001280-77	001260-77	10

SonoBlock II	Size	Item no.	NRFit Item no.	PU
Facet S tip	22G x 50mm (2")	001281-74		10
	22G x 80mm (3 1/5")	001281-71		10
	21G x 100mm (4")	001281-77		10
	20G x 120mm (4 3/4")	001281-72		10

NRFit® is a registered trademark of GEDSA used with their permission.

**PAJUNK® GmbH Medizintechnologie**  
International Sales  
Karl-Hall-Strasse 1  
78187 Geisingen, Germany  
Tel.: +49 7704 9291-0  
Fax: +49 7704 9291-600  
info@pajunk.com  
pajunk.com

**PAJUNK® Medical Produkte GmbH**  
Central Europe  
Pajunkstrasse 2  
78187 Geisingen, Germany  
Tel.: +49 7704 8008-0  
Fax: +49 7704 8008-150  
info@pajunk-medipro.com  
pajunk.de & pajunk.eu

**PAJUNK® Swiss AG**  
Switzerland  
Neuhofstrasse 5A  
6340 Baar, Switzerland  
Tel.: +41 417 609 627  
info@pajunk.ch  
pajunk.ch

**PAJUNK® Medical System L.P.**  
United States of America  
4575 Marconi Drive  
Alpharetta, GA 30005, USA  
Tel.: +1 888 972 5865  
Fax: +1 678 514 3388  
info@pajunk-usa.com  
pajunkusa.com

**PAJUNK® UK Medical Products Ltd**  
United Kingdom  
Unit D1 The Waterfront  
Goldcrest Way, Newburn Riverside  
Newcastle upon Tyne, NE15 8NY, UK  
Tel.: +44 191 264 7333  
info@pajunk.co.uk  
pajunk.co.uk